

Steri Hood

19 Midlands Mills, Cross Hills, Keighley, West Yorkshire BD20 7DT

Tel: (01535) 630776

Fax: (01535) 630827

www.ukhtm01-05.co.uk

sales@ukhtm01-05.co.uk

Don't fall foul of the latest HTM 01 – 05 Sterilisation Rules

**AT LAST AN ECONOMICAL ENVIRONMENT THAT WILL HELP YOU
TOWARDS COMPLIANCE WITH
HTM 01 – 05**



Without recourse to the Expensive and Extensive building works normally associated with full Sterilisation Room Installations

STERI HOOD working for you – helping you with

1. Patient safety
2. Staff safety
3. Helping you ensure that you comply with HTM 01 – 05

STERI HOOD total foot print comes in three different sizes to best suit your application

Steri Hood 2 800mm x 2000mm x 2100mm high

Steri Hood 2.5 800mm x 2500mm x 2100mm high

Steri Hood 3 800mm x 3000mm x 2100mm high

All models consist of – Base Unit

1. Mobile Base frame 40mm x 40mm section mild steel powder coated frame
2. With cabinets of your choice
3. 316 Stainless Steel Rinse sink bowl
4. 316 Acid Resisting Stainless Steel Wash up bowl
5. Waste connections for Disinfector Washer & Autoclave Steriliser

Top Unit

1. Sealed Hood Frame 25mm x 25mm section powder coated frame
2. 6mm back panel and side panels
3. 6mm Perspex sliding doors designed to afford total access to the working base
4. 2no double electrical sockets
5. Sealed fluorescent light to guarantee 800 Lux inside [Steri Hood](#)
6. Dedicated consumer unit
7. Extractor fan designed to continually evacuate the air and any air borne contaminants from within the hood
8. Extract Fan failure audible alarm
9. 100mm Extraction spigot for site connection (optional left - right - back or front)

Options Can Include

1. End mounted wash hand basin
2. Additional under bench cabinets
3. Water Heaters
4. Drainage pumps
5. Mobile used instrument delivery unit

IF YOU MOVE PREMISES YOUR STERIHOO

(YOUR £££ INVESTMENT)

CAN BE EASILY MOVED AND INSTALLED INTO YOUR NEW PREMISES

The whole unit is completely cleanable and easily disinfected

**You may be eligible for grant towards the cost the units. Please
check with you local PGCT**

Steri Hood questions and answers

Q. What do I need to create an environment that will help me to **comply with HTM01-05** and help safe guard my staff and patients?

A. A **Steri Hood** or a dedicated decontamination room

Q. **Where can I install my Steri Hood**

A. **Steri Hood** should be installed out of treatment areas – Away from food preparation areas – & out of general public access areas

Q. **What space will I need to install my new Steri Hood ?**

A. Typical **Steri Hood** 2 Footprint: 2000mm Long x 800mm Deep x 2100mm High

Typical **Steri Hood** 2.5 Footprint : 2500mm Long x 800mm Deep x 2100mm High

Typical **Steri Hood** 3 Footprint : 3000mm Long x 800mm Deep x 2100mm High

Q. **How will I get it into my premises?**

A. **Steri Hood** s can be supplied fully built or in 2No Pieces – Base 900mm x high & Hood 1200mm high In exceptional circumstances where premises offer exceptionally limited access the units can be supplied in component parts and assembled inside your premises.

Q. **What can I have inside my Steri Hood**

A. **Steri Hood** s are built to customers requirements – so for instance your practice might have an autoclave & a disinfecting washer – so you would require 1 dirty wash sink 1 Rinse sink plus the additional plumbing electrical supplies required for both machines – but you may only have only an auto clave – so HTM01-05 recommends that you have a dirty wash sink and a rinse wash sink plus a washer disinfectant & an autoclave – all you need to do is inform the **Steri Hood** sales team of your requirements and we can quote accordingly – if you can provide a rough sketch of your requirement's so much the better – but if not don't worry we will work it out with you.

Q. **I know that also require a hand washing facility - Can I incorporate a Hand washing sink into my Steri Hood**

A. Yes we can install a stainless steel hand washing sink to the right or left hand of the **Steri Hood** dependant on your location – This are positioned outside of the hood - please refer to pictorial illustrations

Q. **Which plumbing services will my Steri Hood Require?**

A. Cold water at minimum

Q. What about hot water?

A. If you have hot water available that ok – if not your **Steri Hood** can be fitted with an instantaneous hot water heater

Q. What about drainage?

A. If you a have localised drain point fine – if not your **Steri Hood** can be fitted with a drain pump

Q. Which electrical services will my **Steri Hood require?**

A. This will depend entirely on your equipment requirement we have typical electrical drawings and calculations these are available on our web site for client down load - **If you are unsure please contact your qualified electrician for advice**

Q. Does my **Steri Hood have to be sited near an outside wall?**

A. No the **Steri Hood** extraction system can be ducted through the building – But internal duct work should always be under negative pressure so the fan must always be sited at the point of duct exit from the building and the external ducting must not discharge under any building opening or into public areas – it is recommended that the extraction ducting should vent above the building eaves .It is always preferable to design your system with the minimum length of ducting possible.

Q. Do you have your own **Steri Hood installation team's?**

A. Yes we have dedicated teams of professional installers - the cost of this service varies and can be quoted on application to mailto:sales@ukHTML01-05.co.uk

Q. Do clients have to use **Steri Hood team of installers?**

A. No - any qualified team of plumbers and electricians should not have a problem installing and maintaining your **Steri Hood**

Q. Do you offer a survey and design service?

A. Yes we do offer a survey and design service – the cost of this service varies and can be quoted on application to mailto:sales@ukHTML01-05.co.uk

Steri Hood

A better way todays compliance requirements

CONTACT STERIHOO NOW

01535 630776

sales@ukhtm01-05.co.uk